



AQUALINE
HEIGHT-ADJUSTABLE BATHROOM EQUIPMENT

For comfort and safety

Coping with the activities in the bathroom may be difficult for elderly and disabled people, because the room is too small and the installations and equipment are fixed. They often sit on a stool or in a wheelchair at the washbasin. This may be difficult if the washbasin is installed to suit a standing user, because of the difference in reach and height. Some elderly and disabled people need help in the bathroom and the fixed installations and equipment will make the work of the helper strenuous because of awkward working postures. Aqualine height-adjustable bathroom equipment will improve the conditions of user as well as helper when working in the bathroom. An Aqualine height-adjustable washbasin will make the user more independent and facilitate the work of the helper. The height-adjustable washbasin offers better working postures for children and adults, seated as well as standing, as the working height may be adjusted to suit the individual user.



40-14600

BASIS WASHBASIN

Basis height-adjustable washbasin includes mirror, shelf, single-lever mixer tap, waste trap, flexible hoses and safety stop under the washbasin.

The height is electrically adjustable from 60 to 100cm by means of a hand control unit with spiral cord. As an alternative to the conventional washbasin, a washbasin for people with disabilities, with or without tilting facility, is available.

40-14600	Basis - 120x64x60cm
40-14620	Basis, with washbasin for people with disabilities - 120x66x75cm
40-14610	Basis, tiltable - 120x66x75cm

ACCESSORIES

40-14690	Light fitting to be mounted above mirror - 60cm
40-14680	Extended lever for mixer tap
40-14685	Non-touch mixer tap
97002452	Large mirror - 70x47cm

(When provided with large mirror or non-touch mixer tap the unit has no shelf)



40-14650

MAXI WASHBASIN

Maxi height-adjustable washbasin includes table, mirror, wall cupboard with shelf, single-lever mixer tap, waste trap, flexible hoses and safety stop under the washbasin.

The cupboard and table may be placed at the right or left-hand side of the washbasin. The height is electrically adjustable from 60 to 100cm by means of a hand control unit with spiral cord.

40-14650	Maxi, washbasin right - 120x94x66cm
40-14660	Maxi, washbasin left - 120x94x66cm

ACCESSORIES

40-14690	Light fitting to be mounted above mirror - 60cm
40-14680	Extended lever for mixer tap
40-14685	Non-touch mixer tap



FLEX-IN WASHBASIN

With the height-adjustable washbasin Flex-In the flexible pipework will always be hidden, resulting in a beautiful and simple finish. The Flex-In washbasin is electrically height-adjustable from 65-95cm. Flex-In with washbasin, tap and flexible plumbing kit all assembled into one unit. Flex-In is available with a standard washbasin or a washbasin for people with disabilities. If preferred, you may fit other bolt mounted washbasins to the height-adjustable unit.

40-14705	Flex-In, with standard washbasin - 57x59cm
40-14706	Flex-In, with washbasin for people with disabilities - 57x72cm
40-14704	Flex-In, height-adjustable unit - 57x14cm

ACCESSORIES

30-67849	Infra-red hand control
40-14680	Extended lever for mixer tap

(Excluding washbasin, tap and plumbing kit)



FLEXI WASHBASIN

Using Flexi Washbasin wheelchair users have unimpeded space up to the wall for better access. The Flexi washbasin is manually or electrically height-adjustable from 68-98cm.

As Flexi is a modular system the manual version may be easily converted into an electrical version to meet changing needs. Elegant covers hiding the height-adjustable legs are available on request in any RAL colour to match or contrast with the wall.

Manual height adjustment

40-14710	Flexi, with washbasin in the middle - 122x51cm
40-14711	Flexi, with washbasin on the right side - 92x52cm
40-14712	Flexi, with washbasin on the left side - 92x52cm

Electrical height adjustment

40-14715	Flexi, with washbasin in the middle - 122x51cm
40-14716	Flexi, with washbasin on the right side - 92x52cm
40-14717	Flexi, with washbasin on the left side - 92x52cm

The products listed above include washbasin, tap and plumbing kit.

ACCESSORIES

30-67849	Infra-red hand control
40-14680	Extended lever for mixer tap
40-14722	Covers
30-67841	Press pad (14x7cm)
30-69160	Safety stop

To eliminate the risk of trapping a person's legs, Ropox always recommends the use of safety stops for electrically operated frames.



AQUALINE BATHTUB

A height-adjustable bathtub is ideal for elderly and disabled people, who need help to take a bath. Getting into and out of the bathtub is facilitated and lifting by means of assistive devices will be easier as the height of the bathtub may be adjusted to suit the lifting equipment. Also the stooped working postures of the helper will be avoided when the working height suits the individual user and helper.

Aqualine bathtub is installed in an electrically height-adjustable frame. The bathtub is provided with mixer tap and hand shower. The height is electrically adjustable from 65 to 95cm by means of a hand control unit with spiral cord. A shower bed with tip-up facility for mounting on the right or left-hand side of the bathtub is available as an optional extra. A bed guard can be mounted for added safety.

Item no.	40-14142	40-14143	40-14140	40-14144	40-14145	40-14146
Width	70cm	70cm	70cm	75cm	80cm	90cm
Length	140cm	150cm	160cm	170cm	180cm	190cm

(Outer measurements)



AQUALINE SHOWER BED

Sometimes elderly and disabled people must be bathed in a lying position. A height-adjustable shower bed facilitates transfer of the person from a bed, chair or stool. The working conditions of the helper will be improved and stooped postures avoided.

Aqualine height-adjustable shower bed is fitted on a wall-mounted bracket and provided with tip-up facility. The height is electrically adjustable from 30 to 100cm by means of a hand control unit with spiral cord. When not used, the shower bed may be tipped up to vertical position to take up less space in the room. The movement is dampened by the gas spring. The shower bed is easy to work with and has a washable lying area. The lying area is 25cm shorter than the shower bed.

A bed guard can be mounted for added safety.

Item no.	40-25520	40-25521	40-25522	40-25523	40-25524	40-25525	40-25526	40-25527	40-25528	40-25529
Width	88cm	88cm	88cm	88cm	88cm	88cm	88cm	88cm	88cm	88cm
Length	100cm	110cm	120cm	130cm	140cm	150cm	160cm	170cm	180cm	210cm

ROPOX

BETTER WAYS TO BETTER DAYS

Ropox A/S

Ringstedgade 221

DK-4700 Naestved

Tel. +45 5575 0500

Fax +45 5575 0550

info@ropox.dk

www.ropox.com

