




Basis Washbasin

Electrical height adjustable washbasin unit
Easy to reach, whether sitting or standing
Ideal for care centres, institutions and private homes

- Range of height adjustment 40cm from 60-100cms
- Adjustment of the washbasin to the required height by means of hand control switch.
- The unit comprises washbasin, mirror, shelf, single-lever mixer tap, trap, flexible hoses and safety stop under the washbasin.
- The safety stop minimises the risk of trapping a person's legs.
- The Basis unit is available with three types of washbasins.
- The basin for people with disabilities has a wide, shallow front for easy access and elbow support which also prevents dripping water.

Basis Washbasin

Basis unit with standard washbasin with mixer tap, overflow, plug, flexible connecting hoses and flexible waste and trap.
Item 40-14600




Standard hand control switch with spiral cord will be provided with the unit
Item 97000527



Flexible water connection and tray provided with the unit
Item 40-25753



Basis unit with washbasin with tilt for people with disabilities. Washbasin without overflow and plug, but with mixer tap, flexible connecting hoses and flexible waste and trap.
Item 40-14610



Accessories


Extended lever for mixer tap
Item 40-14680



Non-touch mixer tap
Item 40-14685



Basis unit with washbasin for persons with reduced mobility. Washbasin without overflow and plug, but with mixer tap, flexible connecting hoses and flexible waste and trap.
Item 40-14620



Large mirror,
70x47cm
Item 40-14695



ROPOX Tips

- With large mirror or non-touch mixer tap the unit has no shelf.
- The safety stop minimises the risk of trapping a person's legs.
- The safety stop must always be accessible.
- The control switch may be placed on the hook provided.
- See that water connection and drain are placed correctly.
- Max. load of frame: 100kg.

ROPOX
BETTER WAYS TO BETTER DAYS

ROPOX A/S
Ringstedgade 221
DK-4700 Naestved
Tel. +45 5575 0500
Fax +45 5575 0550
info@ropox.dk
www.ropox.com