



DOMINO
SLINGS

For comfort and safety

DOMINO SLINGS

MADE OF GREY POLYESTER COMFORT MATERIAL

Domino Slings are produced using a specially manufactured polyester material, which unites comfort with an incredible strength. The material has an emery finish on the inside. This gives the Sling a cotton-like surface, which is both comfortable for the user and at the same time ensures that the Sling does not slide. The material is fire retardant.



DOMINO HYGIENE is incl. safety belt. Suitable for lifting persons from sitting to sitting position. Suitable for hygiene and bathroom purposes.



DOMINO ACTIVE. Supports the body up to shoulder height. Easy to use, as the leg and back sections are extra padded with special foam material. Safety belt as an accessory.



DOMINO COMBI is designed like the Active Sling, but has an extra built-in headsupport, which is hidden in a pocket in the back of the sling. Safety belt as an accessory.



DOMINO CLASSIC. Suitable for persons with little control over their head and body. Suitable for lifting from a lying position. Supports the whole body. Extra foam padding in leg, body and head section – for better comfort.



DOMINO AMPUTATION. The special design of the straps provides optimum support and balance. Safety belt as an accessory.



DOMINO THORAX. Only for stand-up hoist. Unique design and padding. Provides optimum support for stand-up lifting. Can also be supplied without padding under the arms.

DOLPHIN SLING – MADE OF WHITE POLYESTER NET FOR USE IN BATHROOMS

The water can easily pass through the material. Easy to clean.



DOLPHIN HYGIENE is incl. of safety belt. Suitable for hygiene and bathing purposes. Suitable for lifting persons from sitting to sitting position.



DOLPHIN CLASSIC. Suitable for persons with little control over their head and body. Suitable for lifting from a lying position. Supports the whole body.

DOMINO BLUE MADE OF SMOOTH BLUE POLYESTER MATERIAL (NYLON)

The low friction in the material makes it easier to place and remove the sling from the patient.



DOMINO BLUE HYGIENE is incl. of safety belt. Suitable for lifting persons from sitting to sitting position. Suitable for hygiene and bathroom purposes.



DOMINO BLUE ACTIVE. Supports the body up to the shoulder height. Safety belt as an accessory.



DOMINO BLUE CLASSIC. Suitable for persons with little control over their head and body. Suitable for lifting from a lying position. Supports the whole body.

ACCESSORIES



For the **ACTIVE, COMBI** and **AMPUTATION** Slings.

The use of Safety Belts gives the user a feeling of increased safety. They can easily be fixed to the slings.



Detachable 15cm and 30cm long strap loops are available. These can be fitted to the slings if you need increased possibility of individual adjustment and lifting position.

ORDERING INFORMATION

ART NO.	DESCRIPTION
DOMINO SLINGS – FOR STAND-UP HOIST	
25-43011	Domino Thorax Sling with arm-padding. Small, max waist 70cm
25-43012	Domino Thorax Sling with arm-padding. Medium, max waist 100cm
25-43013	Domino Thorax Sling with arm-padding. Large, max waist 130cm
25-43014	Domino Thorax Sling with arm-padding. X Large, max waist 150cm
25-43015	Domino Thorax Sling with arm-padding. XX Large, max waist 170cm
25-48011	Domino Thorax Sling without padding. Small, max waist 70cm
25-48012	Domino Thorax Sling without padding. Medium, max waist 100cm
25-48013	Domino Thorax Sling without padding. Large, max waist 130cm
25-48014	Domino Thorax Sling without padding. X Large, max waist 150cm
DOMINO SLINGS FOR BATH & TOILET	
DOMINO HYGIENE	
25-48021	Domino Hygiene Sling. Small, incl. safety belt, max waist 70cm
25-48022	Domino Hygiene Sling. Medium, incl. safety belt, max waist 110cm
25-48023	Domino Hygiene Sling. Large, incl. safety belt, max waist 125cm
25-48024	Domino Hygiene Sling. X Large, incl. safety belt, max waist 150cm
25-48025	Domino Hygiene Sling. XX Large, incl. safety belt, max waist 170cm
DOMINO BLUE HYGIENE	
25-47011	Domino Blue Hygiene Sling. Small, incl. safety belt, max waist 70cm
25-47012	Domino Blue Hygiene Sling. Medium, incl. safety belt, max waist 110cm
25-47013	Domino Blue Hygiene Sling. Large, incl. safety belt, max waist 125cm
25-47014	Domino Blue Hygiene Sling. X Large, incl. safety belt, max waist 150cm
DOLPHIN HYGIENE	
25-42021	Dolphin Hygiene Sling. Small, incl. safety belt, max waist 70cm
25-42022	Dolphin Hygiene Sling. Medium, incl. safety belt, max waist 110cm
25-42023	Dolphin Hygiene Sling. Large, incl. safety belt, max waist 125cm
25-42024	Dolphin Hygiene Sling. X Large, incl. safety belt, max waist 150cm
DOMINO SLINGS WITH PADDED LEG- AND BACK SUPPORT	
DOMINO ACTIVE	
25-48041	Domino Active Sling. Small
25-48042	Domino Active Sling. Medium
25-48043	Domino Active Sling. Large
DOMINO BLUE ACTIVE	
25-47041	Domino Blue Active Sling. Small
25-47042	Domino Blue Active Sling. Medium
25-47043	Domino Blue Active Sling. Large

ART NO.	DESCRIPTION
DOMINO SLINGS WITH FLEXIBLE HEAD SUPPORT AND PADDED LEG- AND BACK SUPPORT	
DOMINO COMBI	
25-49031	Domino Combi Sling. Small. Optional use of head support
25-49032	Domino Combi Sling. Medium. Optional use of head support
25-49033	Domino Combi Sling. Large. Optional use of head support
DOMINO SLINGS WITH HEAD SUPPORT	
DOMINO CLASSIC	
25-49011	Domino Classic. Small. With padded leg supports
25-49012	Domino Classic. Medium. With padded leg supports
25-49013	Domino Classic. Large. With padded leg supports
25-49014	Domino Classic. X Large. With padded leg supports
DOMINO BLUE CLASSIC	
25-47841	Domino Blue Classic. Small
25-47842	Domino Blue Classic. Medium
25-47843	Domino Blue Classic. Large
DOLPHIN CLASSIC	
25-41011	Dolphin Classic. Small. With padded leg supports
25-41012	Dolphin Classic. Medium. With padded leg supports
25-41013	Dolphin Classic. Large. With padded leg supports
25-41014	Dolphin Classic. X Large. With padded leg supports
AMPUTATION SLINGS WITH PADDED LEG- AND BACK SUPPORT	
DOMINO AMPUTATION	
25-48031	Domino amputations Sling. Small
25-48032	Domino amputations Sling. Medium
25-48033	Domino Amputations Sling. Large
ACCESSORIES FOR DOMINO SLINGS	
25-49111	Safety belt for Domino Active, Combi and Amputations Sling. Small
25-49112	Safety belt for Domino Active, Combi and Amputations Sling. Medium
25-49113	Safety belt for Domino Active, Combi and Amputations Sling. Large
25-49121	Safety belt for Domino Blue Active Sling. Small
25-49122	Safety belt Domino Blue Active Sling. Medium
25-49123	Safety belt for Domino Blue Active Sling. Large
25-40115	Strap loop 15cm. Can be used on all Domino Slings
25-40130	Strap loop 30cm. Can be used on all Domino Slings

ROPOX

BETTER WAYS TO BETTER DAYS

Ropox A/S
 Ringstedgade 221
 DK-4700 Naestved
 Tel. +45 5575 0500
 Fax +45 5575 0550
 info@ropox.dk
 www.ropox.com