




## Maxi Washbasin

Electrical height adjustable washbasin unit  
Easy to reach whether sitting or standing  
Ideal for care centres and private homes

- Range of height adjustment 40cm from 60-100cms.
- Adjustment of the washbasin to the required height by means of hand control switch.
- The unit comprises washbasin with tabletop, mirror, cupboard with shelf, single-lever mixer tap, flexible hoses and safety stop under the washbasin.
- The safety stop minimises the risk of trapping a person's legs.
- The Maxi unit is available with basin to the right or left.

# Maxi Washbasin

Maxi unit with washbasin to the right, with mixer tap, overflow, plug, flexible connecting hoses and flexible waste and trap.  
Item 40-14650




Standard hand control switch with spiral cord will be provided with the unit  
Item 97000527



Flexible water connection and drain will be provided with the unit  
Item 40-25753



Maxi unit with washbasin to the left, with mixer tap, overflow, plug flexible connecting hoses and flexible waste and trap.  
Item 40-14660



## Accessories

Extended lever for mixer tap  
Item 40-14680



Non-touch mixer tap  
Item 40-14685



## ROPOX Tips

- The safety stop minimises the risk of trapping a person's legs.
- The safety stop must always be accessible.
- The control switch may be placed on the hook provided.
- See that water connection and drain are placed correctly.
- Max. load of frame: 100kg

**ROPOX**  
BETTER WAYS TO BETTER DAYS

ROPOX A/S  
Ringstedgade 221  
DK-4700 Naestved  
Tel. +45 5575 0500  
Fax +45 5575 0550  
info@ropox.dk  
www.ropox.com