



Shower bed

Electrical height adjustable shower bed

Folds up to give clear floorspace when not in use


Transfer possible of heights of 30cm-100cm from chair or hoist and reversed

- Adjustment of the shower bed to the required height by means of hand control with cable.
- The shower bed is security fitted to the wall.
- The highest position upto 100cm help carers maintain a good posture.
- The water catchment tray and flexible waste pipe mean less splashing water in the room.
- The shower bed is available in various sizes.


Shower bed

Height adjustable shower beds, range of adjustment 30-100cm.
The net/cover is 25cm shorter than the shower bed.


Item 40-25522




Item 40-25523




Item 40-25524




Item 40-25525




Item 40-25526



Item 40-25527




Item 40-25528



Item 40-25529

Available with cover only



* Lying area

Standard hand control
Will be provided with
the shower bed
Item 97000538



Flexible water connection and trap
Item 40-25752



Accessories

Bed guard 110cm
NOTE: When collapsed
the bed guard will be
longer than the shower
bed for the
lengths 120-130-
140cm
Item 40-25961




Bed guard 60cm,
collapsible
Item 40-25966



Removable shower bed cover with holes and zip.
To be ordered with shower bed.

- Item 40*25520-120
- Item 40*25520-130
- Item 40*25520-140
- Item 40*25520-150
- Item 40*25520-160
- Item 40*25520-170
- Item 40*25520-180



ROPOX
BETTER WAYS TO BETTER DAYS

ROPOX A/S
Ringstedgade 221
DK-4700 Naestved
Tel. +45 5575 0500
Fax +45 5575 0550
info@ropox.dk
www.ropox.com